


Elite Imaging Table



Precision 0047-TM5(H)

*Unobstructed 3-D
c-arm imaging*

500lbs/227kg weight capacity

High Quality Hospital Grade Mobile C-Arm Imaging Table

The "Elite" series of imaging tables are designed for high patient throughput utilizing the highest quality components and designed for ease-of-use by the medical staff with patient comfort in mind.

The tabletop has been designed to give maximum unobstructed imaging allowing 360° rotation radiation-free coverage around the longitudinal axis of the carbon fiber tabletop which makes it ideal for 3D imaging studies.

The tabletop pad is also designed to give the patient maximum comfort being manufactured utilizing our patient form-fitting Visco/Memory "Sofcare" High Density foam.

150mm (58") Full 360° Unencumbered Rotational Imaging Area



Model: 0047-TM3

Multi-purpose, 3-movement
Imaging Table

- Motorized Elevation
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt

Model: 0047-TM4

Extend-Retract, 5-movement
Imaging Table

- Motorized Elevation
- 15.75" (400mm) Motorized Longitudinal Travel
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt





Movement One

- Motorized Elevating Lift
- Lift Travel: (11.75") 300mm
- Optional Travel: (15.75") 400mm



Movement Two

- Motorized Trendelenburg Tilt
- ± 15 Degree Tilt

Movement Three

- Motorized Extend and Retract
- Longitudinal Travel: (15.75") 400mm



Movement Four

- Motorized Lateral Tilt
- Lateral Tilt: ± 15 Degree Tilt



Movement Five

- Motorized Lateral Drive
- Lateral Travel: ($\pm 4'$) ± 100 mm

